



OUTLAB PROCESSING MAIL-IN FORM

For 35mm B&W, 120/220 film
or E-6 slide processing.

For 35mm C41 film please fill out the
"McBain Lab Processing" form available
for download on our website.

LAB USE

REFERENCE

ALL C41 CHECK

TOTAL PUSHES TO INVOICE

DESCRIPTION

Please fill out all fields clearly in print and send it in with your processing order

FULL NAME

EMAIL

ADDRESS

PHONE

CITY/PROVINCE

NOTES

POSTAL CODE

YOUR ORDER:

Number of Rolls	Film Type (35mm/120/220) <small>*Please specify if Xpan/Panoramic 35mm</small>	Scan Size 6MB or 18MB <small>Scanned to DVD unless USB Drive supplied.</small>	Print Size	Finish	Push/Pull <small>Please provide the # of stops. (use + for push, - for pull) Allow 14-21 days for delivery.</small>
			<input type="checkbox"/> 4 x 6 <input type="checkbox"/> 5 x 5 <input type="checkbox"/> 5 x 7	<input type="checkbox"/> Matte <input type="checkbox"/> Glossy	
			<input type="checkbox"/> 4 x 6 <input type="checkbox"/> 5 x 5 <input type="checkbox"/> 5 x 7	<input type="checkbox"/> Matte <input type="checkbox"/> Glossy	
			<input type="checkbox"/> 4 x 6 <input type="checkbox"/> 5 x 5 <input type="checkbox"/> 5 x 7	<input type="checkbox"/> Matte <input type="checkbox"/> Glossy	
			<input type="checkbox"/> 4 x 6 <input type="checkbox"/> 5 x 5 <input type="checkbox"/> 5 x 7	<input type="checkbox"/> Matte <input type="checkbox"/> Glossy	

Most film stock we can pushed up to 3 stops and push/pull up to 1, though this may vary depending on the film stock. We will push or pull as per your instructions, however we advise checking the manufacturer's recommendations regarding push/pull processing for the film. Push processing will increase contrast and grain as well as affect the colour balance.

PRICING - For pricing information please visit <https://www.mcbaincamera.com/film>

SCAN SIZE - Our approximate scan size (shortest edge) is 1200 pixels for 6MB and 2000 pixels for 18MB.

Upon receiving your film we will email you an acknowledgement and invoice, we will also provide an estimated date for delivery of your scans. When ready, a download link will be emailed to you.

MAIL YOUR ORDER TO:

10452 Mayfield Road NW, Edmonton, Alberta T5P 4P4